

MICROPLASTICS FOR BREAKFAST

Mikroplastika za doručak

*Prvi susret istraživača mikro i nanoplastike iz Srbije, Crne Gore & Bosne i
Hercegovine*

Preparation and characterization of microplastics for environmentally relevant laboratory research

Assist. prof. dr. Ula Putar

University of Ljubljana, Faculty of Chemistry and Chemical Technology

Vrčin, 10. 4. 2025



FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

Introduction



University of Ljubljana, Faculty of Chemistry and Chemical Technology, Ljubljana, Slovenia

Chair of Chemical Process, Environmental and Biochemical Engineering



PLANterastic research team:
<https://planterastics.fkkt.uni-lj.si/>



- > 10 years of microplastic research
- > 35 published papers
- 4 book chapters
- 6 national research project
- 1 international research project




Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

Microplastics




Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

Microplastics



20 μm

A microplastic particle...



800 μm

... also a microplastic particle

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>




FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

Microplastics for laboratory research



Preparation



Characterization



Aging/weathering

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



Microplastics for laboratory research



Preparation

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>





FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

1. Microplastics preparation

- Purchase

- Industrial microplastics - specific size and chemical composition
 - Spheres → low environmental relevance
 - Fragments



- Products that contain microplastics (e.g. cosmetic and personal care products) - spherical and fragments, environmentally relevant, characterization needed

<https://youtu.be/RjXNXNHGQys>



Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>




FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

1. Microplastics preparation

- **Preparation**

- Bottom-up (i.e., chemical synthesis) - various shapes, sizes (characterization is needed), environmentally relevant?




Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



1. Microplastics preparation

- Preparation

- Top-down (i.e., fragmentation) - various shapes, sizes (characterization is needed), environmentally relevant



Plastic items
Pellets

Fragmentation

Microplastic
fragments



Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>

Hydrophobic Polymers > 427772

427772 ▶ Sigma-Aldrich

Polyethylene

Write a review Ask a question

average M_n ~4,000 by GPC, average M_w ~1,700 by GPC

Sign In to View Organizational & Contract Pricing

GO Copy Link Email

All Photos (3)

Documents

SDS

COO/COA

Specification Sheet

More Documents (3)

Select a Size

250 G	1 KG
€67.10	€225.00

Change View

Synonym(s): PE

Linear Formula: $\text{HCH}_2\text{CH}_2\text{H}$



PLANT
ERASTICS



FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

1. Microplastics preparation

- Preparation
 - Top-down: Fragmentation



Centrifugal mill



Ball mill



Blender/mill

- Liquid nitrogen
- Dry ice

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

1. Microplastics preparation

- Preparation
 - Top-down: Fragmentation

Example: Centrifugal mill



https://www.youtube.com/watch?v=fR01N7NF4N0&feature=emb_title

Example: Ball mill




<https://www.youtube.com/watch?v=VJDEuF8CEWI>

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



1. Microplastics preparation

- Preparation
 - Top-down: Fragmentation
 - Obtaining particles of specific size? → **Sieving**



<https://www.youtube.com/watch?v=50mmkOS1Fek&list=PLAUyTgal3dz80NvJrpj4DvClpsO2wAxdR&index=7>




Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



Microplastics for laboratory research



Preparation



Characterization




Aging/weathering

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



Microplastics for laboratory research



Characterization

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology





2. Microplastics characterization

- Shape

Light microscope



Stereomicroscope







Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



2. Microplastics characterization

- Shape/Surface morphology

Scanning electron microscopy (SEM)
Transmission electron microscopy (TEM)



Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

2. Microplastics characterization

- **Chemical composition**
 - Fourier-transform infrared spectroscopy (FTIR)
 - Raman spectroscopy
 - Thermal analysis and mass spectrometry (PY-GC-MS)
 - Laser based methods (LIBS, LA-ICP-MS...)
 - Scanning electron microscopy coupled with energy dispersive X-ray spectroscopy (SEM-EDX)
 - ...

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>




FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

2. Microplastics characterization

- Chemical composition

- Fourier-transform infrared spectroscopy (FTIR)



2916 - CH₂ - asymmetric C-H stretching

2848 - CH₂ - symmetric C-H stretching

1464 - CH₂ - bending

719 - CH₂ rocking

- Comparisons with spectral libraries and articles

- FTIR free libraries:

- Spectra Chem:

<https://spectra.chem.ut.ee/>

- Wiley: <https://spectrabase.com/>

- Cameo material database:

https://cameo.mfa.org/wiki/Category:Materials_database

- <https://webbook.nist.gov/chemistry/>

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>




FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology


2. Microplastics characterization

- Chemical composition → elemental composition
 - Scanning electron microscopy coupled with energy dispersive X-ray spectroscopy (SEM-EDX)



PE with adsorbed nTiO₂

PE modified with magnetic fluids



Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>




2. Microplastics characterization

- Size analysis
 - Sieving analysis
 - Optical/stereo microscope
 - Laser diffraction analysis
 - ...



- Mean \pm SD
- Min. - Max.
- **Particle size distribution**




Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



2. Microplastics characterization

- Size analysis
 - Particle size distribution - volume and number






Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



2. Microplastics characterization

- Size analysis
 - Optical microscope



$$5765 \mu\text{m} + 626 \mu\text{m} = 6391 \mu\text{m}$$


Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



2. Microplastics characterization

- Size analysis
 - Photography (digital camera) and processing

Larger particles (>2 mm)



AxioVision program

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>




2. Microplastics characterization

- Size analysis
 - Laser diffraction analyser



<https://youtu.be/ZQqIOrEPI3I>



Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



2. Microplastics characterization





- Are results always accurate?

Small/large particles

Mean \pm SD

Numerical: $23 \pm 7 \mu\text{m}$

Volume: $351 \pm 179 \mu\text{m}^3$




Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>




2. Microplastics characterization

- Are results always accurate?

Agglomeration



Size measurement
with addition of
surfactant



Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

2. Microplastics characterization

- **Other characteristics**
 - Number of particles per mass
 - Specific surface area (BET method)
 - Zeta potential, point of zero charge (DLS - dynamic light scattering, potentiometric titration)
 - Crystallinity (DSC - differential scanning calorimetry, XRD - X-ray diffratography)
 - Wettability (contact angle)
 - Molar mass (chromatography)
 - Density (pycnometer)
 - Presence of additives (GC/LC-MS)
 - ...

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>




FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

Microplastics for laboratory research



Preparation



Characterization



Aging/weathering

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



Microplastics for laboratory research



Aging/weathering

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



PLANT
ERASTICS






FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

3. Microplastics aging

- **Abiotic aging**
 - Accelerated weathering, e.g. Suntest, chemical oxidation
 - Laboratory approach (simulating natural conditions)
 - Outdoor approach



Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>




FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

3. Microplastics aging

- Biotic aging/biofouling




Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



3. Microplastics aging

- **Biotic aging/biofouling**
 - **Natural aging** in the environment (sea, lake...): environmentally relevant but difficult for microplastics (aged in mesh)
 - **In laboratory**
 - Specific monocultures - low environmental relevance
 - **Combined approach:** natural stream water, incubation in laboratory, increased environmental relevance



Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>




FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

3. Microplastics aging

- Biotic aging/biofouling

Combined approach: natural stream water, incubation in laboratory, increased environmental relevance



Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



3. Microplastics

- Biotic agents

Combined application of water, incubation, and increased energy



Microplastics

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>




FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

Conclusions

- Different needs in different fields:
 - Development of analytical methods
 - Ecotoxicology
 - Interaction with other contaminants
 - Fate and behaviour in the environment
 - Degradation
 - ...



But always characterize your particles!

Ula Putar, ula.putar@fkkt.uni-lj.si, <https://planterastics.fkkt.uni-lj.si/>



FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology

Thank you for your attention!



<https://planterastics.fkkt.uni-lj.si/>



<https://si.linkedin.com/showcase/planterastics/>



<https://www.facebook.com/planterastics>



@Planterastics

Dr. Ula Putar

ula.putar@fkkt.uni-lj.si



Dr. Gabriela Kalčíková
Team leader



Dr. Ula Putar
Postdoc



Dr. Tina Skalar
Postdoc



Dr. Gregor Marolt
Postdoc



Barbara Klun
PhD student



Mark Starin
Researcher



Janja Novak
Researcher



Polona Jamnik
Researcher



PLASTouch

plast SENSING))) 



FKKT

UNIVERSITY OF LJUBLJANA
Faculty of Chemistry and Chemical Technology